CL200+

EASY-DIAL PHONE WITH VOLUME AND TONE ADJUSTMENT
TELEPHONE MULTIFONCTIONS A GRANDES TOUCHES

User guide / Notice d'utilisation

E

CAUTION: THIS TELEPHONE PRODUCT IS SPECIFICALLY DESIGNED FOR PEOPLE WITH HEARING DIFFICULTIES AND IS PROVIDED WITH A HANDSET EARPIECE VOLUME CONTROL FOR INDIVIDUAL REQUIREMENTS. DUE CARE MUST BE TAKEN BY ALL USERS THAT THE HANDSET EARPIECE VOLUME CONTROL IS SET TO THE LOWEST LEVEL ACCEPTABLE BY EACH RESPECTIVE USER. CARE SHOULD THEREFORE BE EXERCISED TO ENSURE THAT ANY NAIVE USERS ARE ARE FAMILIAR WITH THIS REQUIREMENT.

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Description	6
Setting up	
- Connecting the telephone	8
- Alarm and ring options	9
- Ring only	
- Ring and shaker	
- Ring and strobe flash	
- Ringer settings	10
- Tone/Pulse switch (T/P)	10
- Time break recall switch (300/100)	11
- Hotel message waiting switch	11
- Wall mounting	12
Using the phone	
- Answering incoming calls	13
- Without headset	
- With headset (option)	
- In handsfree mode	
- Volume and tone receiving	14

- Making a call	16
- Without headset	
- With headset (option)	
- In handsfree mode	
- Last number redial	18
- Recall, * & # buttons	18
- Mute button	18
- Storing telephone numbers	19
- Dialling stored numbers	20

Options:

- Direct Vibe [®] handset	22
- Shaker	23
- Headset with microphone	23

Guarantee	24
-----------	----

For product support and help visit our website at

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DESCRIPTION

Hotel message switch (1)

* Only Geemarc references
can be connected to the
sockets

* Shaker socket (2)
Only ref: SHAKCL_WH

* Headset socket (3)
Only ref: HEADSCL_WH

Tone reception
adjustment (4)

Clip for wall mounting (5)

Speaker volume adjustment (6)

Volume reception
adjustment (7)


Speaker (8)

Last number redial (9)

Time break recall button (10)

Flashing incoming ring indicator (strobe light) (11)

Memory storing key (12)

Tone and volume adjustment automatically activated each time you hang-up the phone (24)

*12 V DC adaptor socket (23)

Line socket (22)

3 priority numbers (21)

Key to activate tone and volume adjustment function (20)

Ringer level adjustment (18)

Ringer select mode (17)

Direct memory buttons (16)


Handsfree / Headset button (15)

Micro

Mute button (13) Hotel message waiting LED (14)

12)

SETTING UP

(*) Classified TNV-3 according to EN60950 standard.

(**) Classified "hazardous voltage" according to EN60950 standard.

1 - CONNECTING THE TELEPHONE

- 1 - Connect the handset cord (31).
- 2 - Open the battery door by unscrewing the retaining screw. Insert 4 x AAA alkaline batteries. Snap the battery door in place.

SETTING UP

In case of power failure, batteries are required for the ring.

- 2 - Connect the line cord (25).
- 3 - Connect the shaker cord (optional).
- 4 - Connect the headset cord (optional).
- 5 - Connect the PTT plug to the wall socket.
- 6 - Connect the 12V DC adaptor cord (27) to the telephone.
- 7 - Connect the 12V DC adaptor (26) to your wall socket (**).

(**) Classified "hazardous voltage" according to EN60950 standard.

2 - ALARM AND RINGER OPTIONS

A 3 position switch (17) located on the right side of the phone allows the settings below:

- 0 : Only the ringer is activated.
-  : Ringer and shaker are activated.
-  : Ringer and strobe flash are activated.

SETTING UP

3 - RINGER SETTINGS

A switch (18) located on right side allows the ringer to be switched HI,LOW,OFF:

 : Off,  : Low,  : Hi

Your phone comes with the ringer switched to «High» and the melody set to «6».

	Low	Mid	High
Melody 1	1	2	3
Melody 2	4	5	6(by default)
Melody 3	7	8	9

If you don't like the ringer melody, you can choose a new one by changing the settings as follow:

- Lift the handset.
- Press  (12) button.
- Press # button.
- Press any button 1 to 9.

Ringer will be activated on next call.

4 - TONE/PULSE SWITCH (T/P)

The Tone/Pulse switch (30) is factory pre-set to

SETTING UP

Tone position (**T**). In the UK, all telephone exchanges now use Tone dialling. If your phone does not dial out, it is probably being used from an older private switch-board (PBX). In this case, the Tone/Pulse switch must be moved to Pulse.

5 - TIME BREAK RECALL SWITCH (300/100)

For the UK, the TBR switch (34) is factory pre-set to 100ms. It can be switched to 300ms for other locations or PBX's.

6 - HOTEL MESSAGE WAITING SWITCH

These are used with the new services provided on hotel network.

A switch  (1) located behind your phone allows to set **ON** the hotel message waiting function.

When the switch is setting **ON** and when you receive a message, the red light  (14) is lit to inform you that you have a message on phone .

SETTINGS

7 - WALL MOUNTING

To wall-mount your telephone, slide out the clip (5), turn around and slide back into the wall mounting position (located in the base unit, just below the ear piece). This will ensure that the handset stays firmly in position. Make 2 holes in the wall 80mm apart (use the wall mount template) , insert the wall plugs and the wall mounting screws. The screws should protrude from the wall by 6-7 mm. Place the phone onto the screw-heads (32) and slide down to secure.

If you are using this telephone with a hearing aid, set your hearing aid to the “T” setting.


USING THE PHONE

1 - ANSWERING INCOMING CALL



Without headset:

- When an incoming call is received, the phone will ring, to answer the call lift the handset and speak.
- On completion of the call, carefully replace the handset in the cradle to release the line.

With headset (option):

- When an incoming call is received, and if you wish to answer the call with the headset, press the  (15) button and speak.
- On completion of the call, press the  button again to release the line.

In handsfree mode:


- When an incoming call is received, the phone will ring, to answer the call press the  button.
- On completion of the call, press the  button to release the line.

USING THE PHONE


2-VOLUME AND TONE RECEIVING ADJUSTMENT


You can increase the volume and adjust the receiving tone according to your own hearing requirements. The button →) (20) allows extra-receiving and tone functions to be activated or not during the conversation.

When the function is activated, the red light of extra receiving and tone function located near the button →) is lit. Adjust the **TONE** (4) and **VOLUME** (7) buttons to remove the level according to your own hearing requirements.

USING THE PHONE

Note: When the tone (treble), earpiece volume and sending levels are set to the maximum, whistling may occur. If whistling occurs the sending level settings should be decreased (switch  (19)). When the switch  (24) located behind the phone is **ON**, the extra-receiving and tone functions are automatically activated each time you use the phone. When the switch  is switched **OFF**, you will need to activate the function by pressing the button  (20). Tone and volume adjustment (15dB to 30dB) is only activated when the button is set to **ON**. When this button is set to **OFF** only the volume can be adjusted from 0dB to 15dB.


USING THE PHONE

3 - MAKING A CALL


Without headset:

- Lift the handset.
- Await the dial tone and dial required number.
- On completion of the call, carefully replace the handset in the cradle.

With headset (option):

- Press the  (15) button.
- Await the dial tone and dial required number.
- On completion of the call, press the  button again to clear the line.

In handsfree mode

- Press the  (15) button.
- Await the dial tone and dial required number.
- On completion of the call, press the  button to release the line.
- If your correspondent doesn't answer your call, press the  button again to release the line.

USING THE PHONE

Note:

- To deactivate the handsfree function during a call, lift the handset.

To activate the handsfree function during a call, press the  button and place the handset on the cradle.

- During a call in handsfree, you can adjust the speaker volume to a desirable level.

Simply use the slide control VOLUME (6) located on the left of the unit.

USING THE PHONE

4 - LAST NUMBER REDIAL

- Lift the handset.
- Await the dial tone and press  (9) button.
- The previously dialled number will be automatically redialled (not in use for memory numbers).

5 - RECALL, * & # BUTTONS

These are used with the new services provided on digital exchanges. For details please contact your network operator.

6 - MUTE BUTTON

If you wish to speak privately to someone else in the room without your caller overhearing, press and hold the  button (13). You will still be able to hear the callers but they cannot hear you during the mute operation. To resume normal conversation, release the  button.

USING THE PHONE

7 - STORING TELEPHONE NUMBER

Your phone can store 12 direct memories including 3 priority numbers (21 digits maximum for each memory).

For example, the 3 priority **P1**, **P2**, **P3** memories (21) allow you to store the 3 numbers below:

- Local, Police station,
- Your doctor
- 1471

It is recommended to store these numbers during installation

To store a number please follow instructions below:

- Lift the handset.
- Press  button (12).
- Press **M1**, **M2**, ..., or **M9** (16) or one of 3 priority numbers (21) to store your number.
- Dial the number.
- Press  button again.

Numbers will be lost after a few hours if you disconnect your phone. Entering a new number in memory will automatically erase the previous number.

USING THE PHONE

8 - DIALLING STORED NUMBER

- Lift the handset and await the dial tone.
- Press appropriate button **M1**, **M2**, ..., or **M9** (16) or one of the 3 priorities **P1**, **P2**, **P3** numbers (21), the stored number will be dialled automatically.

To register your stored phone numbers, use the index card located under the base.

Use this table below to prepare the phone book:

	1	
	2	
	3	
	4	
	5	
	6	
	7	
	8	
	9	

Example :

Bruce	1	01 64 57 ...
-------	---	--------------

USING THE PHONE

For product support and help visit our website at

www.geemarc.com

telephone 01707 384438

or fax 01707 372529

OPTIONS

1- DIRECT VIBE ® HANDSET (OPTIONAL)

If you have conductive hearing loss you can use a Direct Vibe Handset, try placing the Direct Vibe handset on your temple, in front and slightly below your ear or behind your ear to see which position works best for you.

If you don't have hearing loss you can use your phone normally and put the receiver on your ear.

If you require a Direct Vibe Handset ref: HSDVCL_WH you may purchase it from your retailer or audiologist.

In case of difficulty please visit our website (www.geemarc.com) or phone our HELPLINE 01707 384 438 .

To set up the Direct Vibe handset only disconnect the standard handset on the handset socket (31) then replace it with the Direct Vibe Handset.

OPTIONS

2- SHAKER (OPTIONAL)

If you wish to complement you phone by using a **SHAKER** ref: SHAKCL_WH connect the shaker plug in the socket  (2) located behind your phone. When an incoming call is received, the shaker will be activated during ringing.

Note: This option must be set in the ALARM AND RINGER OPTIONS item (page 9) to be activated. To purchase a SHAKER see item 1 above.

3- HEADSET WITH MICROPHONE (OPTIONAL)

If you wish to answer or make calls by using a Geemarc headset ref: HEADSCL_BLK, connect the headset plug in the headset socket  (3) located behind your phone.

Note: When using a headset, the speaker will be activated and it's recommended to adjust the speaker (6) to the low position.

GUARANTEE

GUARANTEE

From the moment your Geemarc product is purchased, Geemarc guarantee it for the period of one year. During this time, all repairs or replacements (at our option) are free of charge. Should you experience a problem then contact our helpline or visit our website at www.geemarc.com. The guarantee does not cover accidents, negligence or breakages to any parts. The product must not be tampered with or taken apart by anyone who is not an authorised Geemarc representative. The Geemarc guarantee in no way limits your legal rights.

IMPORTANT: YOUR RECEIPT IS PART OF YOUR GUARANTEE AND MUST BE RETAINED AND PRODUCED IN THE EVENT OF A WARRANTY CLAIM

Please note: The guarantee applies to the United Kingdom only.

DECLARATION: Hereby Geemarc Telecom SA declares that this product is in compliance with the essential requirements and other relevant provisions of the Radio and Telecommunications Terminal Equipment Directive 1999/5/EEC and in particular article 3 section 1a, 1b and section 3. The telephone does not operate if the line current is lower than 18 mA.

Electrical connection : The apparatus is designed to operate from a 230V 50Hz supply only. (Classified as «hazardous voltage» according to EN60950 standard).

The apparatus does not incorporate an integral power on/off switch. To disconnect the power, either switch off supply at the mains power socket or unplug the AC adaptor. When installing the apparatus, ensure that the mains power socket is readily accessible.

Telephone connection : Voltages present on the telecommunication network are classified TNV-3 (Telecommunication Network Voltage) according to the EN60950 standard.

Recycling directives

The WEEE (Waste Electrical and Electronic Equipment) has been put in place for the products at the end of their useful life are recycled in the best way.

When this product is finished with, please do not put it in your domestic waste bin.

Please use one of the following disposal options:

- Remove the batteries and deposit them in an appropriate WEEE skip. Deposit the product in an appropriate WEEE skip.
- Or, hand the old product to the retailer. If you purchase a new one, they should accept it.

Thus if you respect these instructions you ensure human health and environmental protection.




For product support and help visit our website at

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